

Q-Fit Back Armour Brace

E-WAR021 Q-fit back armour

E-WAR021A Premium Q-fit back armour

Q-fit Back Armour forms to the patient's body for the best back pain treatment and comfort. We utilize the Q-fit system for stability and support to the L1-L5 or T9-S1 lumbar back region. It is also designed to provide a quick and easily-controlled compression with its single dial design. Compared to many back brace, the Q-fit back armour is thinner, more lightweight and breathable. The enhanced stability is backed up by the posterior and side panels that increase support and help alleviate pain and discomfort.

Indications

chronic or acute low back pain, lower back sprain/ strain, disc herniation, degenerative changes, post- op support.

Size (cm)

S	M	L	XL	XXL
56-71	69-84	81-96	94-109	107-122

Optional

Posterior panel	E-WAR021-PP	S/M L/XL
Anterior panel	E-WAR021-AP	ONE SIZE
Side panel (a pair)	E-WAR021-SP	ONE SIZE



E-WAR021A
Premium Q-fit back armour



E-WAR021
Q-fit back armour

Features

- The Q-fit system allows you to do quick, one-handed adjustment so you can control the level of support and comfort for standing or sitting
- 3:1 ratio, low-friction Q-fit System provides powerful, smooth and easily-controlled compression and support.
- Single dial front closure provides easy loosening and tightening
- The low profile posterior panel, anterior panel and side panels all are connected via a Velcro closure, and can be removed and reattached easily.
- Thin, lightweight, breathable mesh material combines with contouring panels for increased support.




The contoured posterior panel with comfortable foam provides superior sagittal support from T9-S1 and great fit.



The contoured anterior panel with comfortable foam provides superior sagittal support from L1-L5 and great fit.



The side panels made of comfortable foam, provide additional lateral control.



The Q-fit system provides an efficient, powerful compression and support for lumbar region. Single dial makes quicker and easier using.

-supported by **Boa Closure System**



elife

Q-Fit Back Armour Brace

E-WAR021

Q-fit back armour

E-WAR021A

Premium Q-fit back armour

e-life International Co., Ltd

7F, No. 1, Baosheng Rd.,
Yonghe Dist., New Taipei City 234,
Taiwan (R.O.C)
TEL: +886-2-2923-1516
FAX: +886-2-2929-0783
EMAIL: info@e-lifebracing.com

embrace life.



Armour Cervical Collar

E-CO007

Armour Cervical Collar forms to the patient's neck for the best neck pain support and care. Compared to many cervical brace, the Cervical Armour simplifies the complexity during applying the brace for saving time and simply using . It also reduces inventory and saves money by one size design. Armour Cervical Collar utilizes unique design and mechanism to enhance patient care efficiently and economically. It will be a nice solution to accommodate patient's cervical condition.

Features


- *Two-piece polyethylene design enhances patient comfort and is easily adjusted with buckle design
- *The unique adjustment technology allows 7 levels of chin height adjustment for customized fitting with a simple knob sliding.
- *Designed to help restrict cervical flexion, extension, lateral bending, and rotation
- *Removable liner is made of soft and moisture wicking coolmax® for greater comfort and proper hygiene.
- *The large tracheal opening provides functional access for the airway management and other essential procedures
- *Safe for use with X-ray,CT, and MRI
- *One size fits all

Indications

- *Post-op immobilization
- *Degenerative disorder
- *Cervical spinal stenosis
- *Cervical sprain and strains
- *Whiplash syndrome
- *Halo removal

Size

One size fits all
(neck circumference : 30-45 cm
Height :11-17.5cm)



Lining Coolmax® wicks moisture away from the body to ensure proper patient hygiene and comfort.



Buckle design makes easy application and saves more time. Only need to adjust desired size at the first application by Velcro tabs



Back panel with ventilation holes improve breathability



The unique adjustment technology allows 7 levels of chin height adjustment for customized fitting with a simple knob sliding



elife

E-CO007

Armour Cervical Collar

embrace life.